

Assemble Lab Cylindrical Press Mold For Laboratory Use

Item Number: PMAC

Introduction

Premium lab cylindrical press mold for flawless sample prep. Prevents delamination, ultra-durable Japanese steel. Custom sizes available. Get yours now!

[Learn More](#)



Model	PMAC-A	PMAC-B	PMAC-C	PMAC-D	PMAC-E
Material	Cr12MoV				
Indenter hardness	HRC60-HRC62				
Sample size	Φ3~Φ4~Φ5~Φ6~Φ8~Φ10mm (M)	Φ12~Φ13~Φ15~Φ18~Φ20mm (M)	Φ30~Φ40mm (M)	Φ50~Φ60mm (M)	Φ70~Φ80 mm (M)
Cavity depth	30mm (N)	40 mm (N)	50 mm (N)	55 mm (N)	60 mm (N)
External dimensions	Φ43*93mm(L*H)	Φ53*120mm(L*H)	Φ73*133mm~Φ95*133mm(L*H)	Φ115*150mm~Φ127*150mm (L*H)	Φ153*180mm~Φ180*180mm(L*H)
Weight(Kg)	0.75Kg	1.2Kg	3.8Kg~6.3Kg	14Kg~20Kg	30Kg~40Kg

Other sizes can be customized