

Assemble Square Lab Press Mold For Laboratory Use

Item Number: PMAS



Introduction

KINTEK's Assemble Lab Press Mold ensures precise sample prep for delicate materials, preventing damage with quick-disassembly design. Ideal for thin strips & reliable demolding.

[Learn More](#)

Model	PMAS-A	PMAS-B	PMAS-C	PMAS-D	PMAS-E
Material	Cr12MoV				
Indenter hardness	HRC60-HRC62				
Sample size	3*3□4*4□5*5□6*6□8*8□10*10 mm (M)	12*12□15*15□18*18□20*20mm(M)	30*30□40*40 mm (M)	50*50□60*60 mm(M)	70*70□80*80 mm (M)
Cavity depth	30mm (N)	40 mm (N)	50 mm (N)	55 mm (N)	60 mm (N)
External dimensions	Φ53*120mm(L*H)	Φ73*133mm(L*H)	Φ95*133mm□Φ115*133mm(L*H)	Φ127*150mm□Φ153*150mm (L*H)	Φ180*180mm□Φ200*180mm(L*H)
Weight(Kg)	1.2Kg	3.6Kg	7Kg□14Kg	20Kg□30Kg	40Kg□50Kg

Other sizes can be customized