

Carbide Lab Press Mold For Laboratory Sample Preparation

Item Number: PMW



Introduction

Premium carbide lab press molds for precise sample preparation. Durable, high-hardness YT15 material, customizable sizes. Ideal for XRF, battery research & more.

[Learn More](#)

Model	PMW-A	PMW-B	PMW-C				
Material	Carbide YT15						
Indenter hardness	HRC85-HRC90						
Sample size	φ3 /φ4 /φ5 /φ6 /φ8 /φ10 mm	φ12 /φ13 /φ15 /φ18 /φ20 mm	φ22 /φ25 /φ28 /φ30 mm				
Cavity depth	30 mm	40 mm	45 mm				
External dimensions	φ43×93 mm	φ53×120 mm	φ73×133 mm				
Weight(Kg)	0.78	1.8	3.8				

Other sizes can be customized