

Heated Hydraulic Press Machine With Heated Plates For Vacuum Box Laboratory Hot Press

Item Number: PCVM



Introduction

Enhance lab precision with KINTEK's Heated Vacuum Lab Press for uniform, oxidation-free samples. Ideal for sensitive materials. Get expert advice now!

[Learn More](#)

Instrument model	PCVM-10T
Pressure Range	0-10.0 tons
Pressurization process	Manual pressurization
Heating temperature	Maximum 500°C
Heating plate	180x180mm, 200x200mm (customizable options may be available)
Vacuum box material	SUS 304 (stainless steel)
Studio Size	400x400x400mm
Door size	300x350mm
Vacuum degree	-0.1MPa
Dimensions (LxWxH)	450x550x850mm
Power supply	220V 50Hz (110V option available)