


Lab Anti-Cracking Press Mold

Item Number: PML



Introduction

Precision Anti-Cracking Press Mold for lab use. Durable Cr12MoV steel, high-pressure resistant, customizable sizes. Ideal for material testing. Get yours now!

[Learn More](#)

Instrument model	PMT
Sample shape	
Die Material	Alloy tool steel: Cr12MoV
Indenter hardness	HRC60-HRC62
Sample size (Standard Options)	Φ6mm, Φ8mm, Φ10mm, Φ15mm, Φ20mm (Customizable - M)
Cavity depth	40mm (Customizable - N)
Dimensions (Outer)	Φ98*120mm (L*H)
Weight	Approx. 5Kg
Diagram of hydraulic powder press size compatibility	