

Lab Cylindrical Press Mold With Scale

Item Number: PCMC



Introduction

KINTEK's Cylindrical Press Mold ensures precision material processing with uniform pressure, versatile shapes, and optional heating. Ideal for labs and industries. Get expert advice now!

[Learn More](#)

Instrument model	PCMC
Sample shape	
Die heating	Alloy tool steel :Cr12MoV
Indenter hardness	HRC60-HRC62
Sample size	Φ10Φ12Φ13Φ15Φ18Φ20 mm (M)
Cavity depth	100mm (N)
Dimensions	Φ53*220mm(L*H)
Weight	4.8Kg
Diagram of hydraulic powder press size	