


Lab Heat Press Special Mold

Item Number: PCHF



Introduction

Precision KINTEK lab press molds for reliable sample prep. Durable, customizable, and ideal for diverse research needs. Enhance your lab's efficiency today!

[Learn More](#)

Instrument model	PCHF
Sample shape	
Die heating	0°C-500°C
Indenter hardness	SUS 304
Sample size	Rectangular or bone shape
Cavity depth	0.75mm□1.35mm□1.75mm□2.75mm
Dimensions	80x80□180x180□200x200mm
Weight	0.4kg□0.8kg□1.0kg
Diagram of hydraulic powder press size	