

Lab Infrared Press Mold For Laboratory Applications

Item Number: PMID

Introduction

KINTEK's lab press molds ensure precise sample preparation with durable tungsten carbide construction. Ideal for FTIR, XRF, and battery research. Custom sizes available.

[Learn More](#)



Instrument model	PMID
Sample shape	
Die material	Tungsten carbide
Indenter hardness	HRC68-HRC85
Standard Sample size	Φ13 mm(M)
Cavity depth	20mm(N)
Dimensions	Φ43*78mm(L*H)
Weight	0.76Kg
Diagram of hydraulic powder press die size	