

Lab Infrared Press Mold For No Demolding

Item Number: PMI



Introduction

Streamline infrared sample prep with KINTEK's non-demountable molds—achieve high transmittance without demolding. Ideal for spectroscopy.

[Learn More](#)

Instrument model	PMI-A	PMI-B
Sample shape		
Die material	Tungsten carbide	
Indenter hardness	HRC68-HRC85	
Sample size	Φ13 mm(M)	Φ7 mm(M)
Cavity depth	10mm(N)	5mm(N)
Dimensions	Φ76*50*70mm(L*W*H)	Φ76*30*70mm(L*W*H)
Weight	0.76Kg	0.35Kg
Diagram of hydraulic powder press size		