

Lab Manual Microtome Slicer For Tissue Sectioning

Item Number: PT10



Introduction

Precision Manual Slicer for labs: Achieve accurate, consistent sample preparation with adjustable high-precision cutting. Ideal for research, food, and industry.

[Learn More](#)

Instrument model	PT10
Mold material	3Cr13 (quenched stainless steel)
Slice size	Φ10, Φ12, Φ14, Φ15, Φ16, Φ18, Φ19(d)
Overall dimensions	120*200*280mm (L*W*H)
Equipment weight	10kg
Machine size diagram	