

Lab Ring Press Mold For Sample Preparation

Item Number: PMO



Introduction

High-precision Ring Press Dies for uniform pellets in labs & industry. Durable Cr12MoV alloy, sizes $\Phi 3$ -80mm. Boost efficiency & accuracy today!

[Learn More](#)

Instrument model	PMQ	
Sample shape		
Die Material	Alloy tool steel: Cr12MoV	
Indenter hardness	HRC60-HRC62	
Sample size options	$\Phi 7$ -3, $\Phi 10$ -5, $\Phi 20$ -10 mm (M)	$\Phi 30$ -10, $\Phi 50$ -20mm (D)
Cavity depth	40mm (N)	45mm (D)
Dimensions	$\Phi 53 \times 120$ mm (L*H)	$\Phi 72 \times 100$ mm, $\Phi 88 \times 120$ mm (D*L)
Weight	1.4Kg	3.5kg, 5kg
Diagram of hydraulic powder press size		