

Lab Round Bidirectional Press Mold

Item Number: PMSY



Introduction

Precision Round Bidirectional Press Mold for lab use, high-density compaction, Cr12MoV alloy steel. Ideal for powder metallurgy & ceramics.

[Learn More](#)

Instrument model	PMSY
Sample shape	
Die material	Alloy tool steel : Cr12MoV
Indenter hardness	HRC60-HRC62
Sample size	Φ12□Φ13□Φ15□Φ18□Φ20mm(M)
Cavity depth	40mm (N)
Dimensions	Φ88*175mm(L*H)
Weight	3.0Kg
Diagram of hydraulic powder press size	