

Laboratory Hydraulic Press 2T Lab Pellet Press For Kbr Ftir

Item Number: KT-KBP



Introduction

KINTEK 2T Lab Hydraulic Press for precise FTIR sample prep, durable KBr pellet creation, and versatile material testing. Ideal for research labs.

[Learn More](#)

Model	PCKBR-2T
Working Pressure (T)	0-2 (30Mpa)
Piston Diameter	Φ45mm
Pressure Stability	≤1MPa/10min
Workbench Diameter	Φ45mm
Number of Columns	Two
Working Space (W×D)	54×55 mm
Dimensions (W×L×H)	100×220×220 mm
Weight (kg)	4.8 kg