

# Special Shape Lab Press Mold For Laboratory Applications

Item Number: PMT



## Introduction

Special Shape Press Molds for precise lab applications. Customizable, high-pressure performance, and versatile shapes. Ideal for ceramics, pharmaceuticals, and more. Contact KINTEK today!

[Learn More](#)

Instrument model	PMT
Sample shape	
Die material	Alloy tool steel: Cr12MoV
Indenter hardness	HRC60-HRC62
Sample size (diameter)	Φ6, Φ8, Φ10, Φ15, Φ20mm (M) (Customizable)
Cavity depth	40mm (N) (Customizable)
Dimensions (Outer)	Φ53*120mm (L*H) (Varies by customization)
Weight	Approx. 1.4Kg (Varies by customization)
Diagram of hydraulic powder press size compatibility	