

Square Bidirectional Pressure Mold For Lab

Item Number: PMS-F



Introduction

Achieve high-precision powder molding with KINTEK's Square Bidirectional Pressure Mold for superior lab results. Explore now!

[Learn More](#)

Instrument model	PMSY
Sample shape capabilities	
Die material	Alloy tool steel : Cr12MoV
Indenter hardness	HRC60-HRC62
Standard Sample sizes (Customizable)	12*12, 15*15, 18*18, 20*20 mm (M)
Cavity depth	40mm (N)
Dimensions	Φ88*175mm(L*H)
Weight	3.0Kg
Diagram of hydraulic powder press die components	