

Xrf Kbr Steel Ring Lab Powder Pellet Pressing Mold For Ftir

Item Number: PMXS

Introduction

Precision XRF steel ring pellet mold for lab sample prep. Durable, efficient, ensures accurate XRF analysis. Custom sizes available. Order now!

[Learn More](#)



Instrument model	PMXS
Sample shape	
Die material	Alloy tool steel: Cr12MoV
Indenter hardness	HRC60-HRC62
Sample size	Φ32mm, Φ40mm (M)
Cavity depth	45mm (N)
Dimensions	Φ73*133mm (L*H)
Weight	3.2Kg
Diagram of hydraulic powder press size	